

The ISO 9001:2008 international quality certification system is adopted by our company

TO BE A WORLD-LEADING ANALYTICAL TESTING SOLUTIONS PROVIDER

- Spectroscopy
- Chromatography
- Mass Spectrometry



EDX1800B X-ray Fluorescence Spectrometer

RoHS Testing | Full-element Analysis for minerals | Plating Thickness Testing | Precious Metals Testing

United Kingdom:

Optech Solutions Ltd.

Riverside Court, Beaufort Park, Chepstow NP16 5UH, UK

Tel: +44 (0)1291 418148, Fax: +44 (0)1291 418143

Website: www.optechsolutions.co.uk

E-mail: info@optechsolutions.co.uk

Test data in this manual, if not noted, is our company's test data.
All information in this manual is for reference only, which is subject to any change without notice.

Skyray Instrument Copyright 2010
Press date:2011.07.05

Instrument Introduction:

With the widespread of EDX1800B over different fields, we design this type to meet the need of optimizing the product performance and improve the safety protection grade.

The reliability of the product is improved by using the high voltage source and X-ray tube of the new generation and the testing efficiency is improved by the adopting the high power of X-ray tube.



Testing 75 kinds of elements, 1ppm limit of detection, Repeatability 0.1%, Stability 0.1%



RoHS



Application fields:

- RoHS testing
- Mining and alloy (Cu, stainless steel and so on) componential analysis
- Measurement of plating thickness, measurement of electroplate liquid and plating content
- The content test of precious metal such as gold, platinum and silver and different kinds of jewelry
- Mainly applied in RoHS directive industries, precious metals and jewelries processing industries; banks, jewelry shops and test institutes; electroplating industries

Good shielding action of X-ray tube of new generation, radiation level of X-ray is equal to that of common atmospheric environment

The performance is stable and reliable, achieving higher test efficiency

The automatic function of door sensor and high voltage lock gives you protection from all directions

Performance advantages:

- Down-side X-ray Source: meet the test requirements of samples of different kinds and shapes
- Collimator and filter: the Auto-switch between various collimators and filters to meet the application of different testing methods
- Movable platform: sophisticated manual movable platform is convenient for locating test point
- High-resolution detector: improve the analyzing accuracy
- High voltage source and X-ray tube of the new generation: the performance is stable and reliable, achieving higher test efficiency

Technical specifications:

- Measurable elements: S to U.
- Limit of detection (LOD) reaches 1ppm.
- Element content: ppm to 99.99%
- Arbitrary optional analysis and identification models
- Independent matrix effect correction models
- Multi-variable non-linear regression procedure
- Excellent repeatability: 0.1%
- Long-time working stability: 0.1%
- Ambient temperature: 15°C to 30°C
- Power supply: AC 220V±5V, AC purified stabilized voltage power supply.
- Energy resolution: 160±5eV
- Sample chamber size: 439mm×300mm×50mm
- Instrument size: 550mm×410mm×320mm
- Instrument weight: 45kg

Standard configurations:

- Movable sample platform
- Signal-to-Noise Enhancer (SNE)
- Electric-cooling Si-PIN detector
- Signal detection electronic circuit
- High and low voltage power
- X-ray tube of high power
- Computer and ink-jet printer

Streamlined man-machine design, promising your operation security


Operation indicator makes you operate comfortably



Test cases:




RoHS Testing




Test Results Spectrum




Mining Testing



Plating Thickness Testing



Precious Metals Testing



What are RoHS and WEEE Directives?

The European Union has adopted Directive 2002/95/EC on the restriction of certain hazardous substances (RoHS) and Directive 2002/95/EC on waste electrical and electronic equipment (WEEE) with their publication in the Official Journal of the European Union on February 13, 2003. WEEE comes into effect on August 13, 2005 and RoHS requires the substitution of various heavy metals (lead, mercury, cadmium and hexavalent chromium) and brominated flame retardants (polybrominated biphenyls [PBB] or polybrominated diphenyl ethers [PBDE]) in new electrical and electronic equipment put on the market from July 1, 2006.



Testing standard of substances restricted by RoHS Directive

Hazardous substances	Standards (mg/kg)
Cd	100
Pb	1000
Hg	1000
Cr⁶⁺	1000
PBB_s	1000
PBDE_s	1000

Restricted substances and their typical uses

Pb	
Solders	
Paints	Pigments and driers
Glass materials	Pb is allowed in fluorescent lamp
Ceramic materials	Pb is allowed in certain electronic ceramic materials
Iron, aluminum and copper materials	A certain amount of Pb is allowed
Plastics	PVC stabilizer and pigments
Batteries	Pb is allowed in acidic batteries for vehicles
Cd	
Plastics	Stabilizer and pigments
Solders	Seldom used
Ceramics	Seldom used
Connectors	Relays and switches
Batteries	Cd is allowed in Ni-Cd batteries
Semiconductors	Optical sensors and solar cell panels
Hg	
Batteries	Prohibited (see battery directive)
Connectors	Relays and sensitive switches
Fluorescent lamps	A certain amount of Hg is allowed
Cr ⁶⁺	
Passivation layers	Commonly used for naked metal surfaces to enhance adhesion of plating layers
Anti-corrosive plating layers	Painting and plating layers
Chrome plating layers	Plating layer of chromium metal is not under control
Plasticizer	Commonly used to plastics plating process but not final products
PBB _s & PBDE _s	
Plastics	Brominated flame retardants

The analytical method of filtration for XRF to test RoHS substances



Characteristic X-radiation of element

Each element will emit X-ray at its own energy level when excited. This X-ray is characteristic and called X-ray fluorescence. It is the foundation of analysis.

Scattering

It is the background of spectrum.

Photoelement

The photoelectron is the foundation of detector. In the sample, the X-ray intensity of every element is expressed as $I_1, I_2, I_3, I_4, I_5 \dots$ respectively. The element content C is the function of X-ray fluorescence intensity I , expressed as follows:

$$C = f(I_1, I_2, I_3, I_4, I_5, \dots)$$

This equation is too complicated and can be simplified as:

$$C = K_1 I_1 + K_2 I_2 + K_3 I_3 + K_4 I_4 + K_5 I_5 \dots$$

Where

C is the element content in the sample; $I_1, I_2, I_3, I_4, I_5 \dots$ are X-ray intensity of element respectively; $K_1, K_2, K_3, K_4, K_5 \dots$ are coefficients which can be determined by measuring known standard sample to calibrate.

PERIODIC TABLE OF ELEMENTS

1	IA 1 H 氢 Hydrogen 1.008	IIA 2 He 氦 Helium 4.008																	0																																																																																				
2	3 Li 锂 Lithium 6.94 0.052	4 Be 铍 Beryllium 9.012 0.110	5 B 硼 Boron 10.81 0.185	6 C 碳 Carbon 12.01 0.282	7 N 氮 Nitrogen 14.01 0.392	8 O 氧 Oxygen 15.99 0.523	9 F 氟 Fluorine 18.99 0.677	10 Ne 氖 Neon 20.17 0.861													18 Ar 氩 Argon 39.94 2.857	19 K 钾 Potassium 39.1 3.51	20 Ca 钙 Calcium 40.08 3.589	21 Sc 钪 Scandium 44.96 3.31	22 Ti 钛 Titanium 47.90 4.088	23 V 钒 Vanadium 50.94 4.459	24 Cr 铬 Chromium 51.99 4.949	25 Mn 锰 Manganese 54.94 5.411	26 Fe 铁 Iron 55.84 5.895	27 Co 钴 Cobalt 58.93 6.400	28 Ni 镍 Nickel 58.7 6.825	29 Cu 铜 Copper 63.54 7.472	30 Zn 锌 Zinc 65.38 8.041	31 Ga 镓 Gallium 69.72 8.631	32 Ge 锗 Germanium 72.5 7.472	33 As 砷 Arsenic 74.92 8.041	34 Se 硒 Selenium 78.9 8.631	35 Br 溴 Bromine 79.90 10.36	36 Kr 氪 Krypton 83.8 11.91	37 Rb 铷 Rubidium 85.47 13.38	38 Sr 锶 Strontium 87.62 14.97	39 Y 钇 Yttrium 88.91 16.75	40 Zr 锆 Zirconium 91.22 17.69	41 Nb 铌 Niobium 92.91 18.65	42 Mo 钼 Molybdenum 95.94 19.63	43 Tc 锝 Technetium (99) 20.17	44 Ru 钌 Ruthenium 101.0 21.69	45 Rh 铑 Rhodium 102.9 22.76	46 Pd 钯 Palladium 106.4 23.86	47 Ag 银 Silver 107.9 24.99	48 Cd 镉 Cadmium 112.4 26.14	49 In 铟 Indium 114.8 27.38	50 Sn 锡 Tin 118.6 28.60	51 Sb 锑 Antimony 121.7 29.85	52 Te 碲 Tellurium 127.6 31.13	53 I 碘 Iodine 126.9 33.78	54 Xe 氙 Xenon 131.3 35.46	55 Cs 铯 Cesium 132.9 36.55	56 Ba 钡 Barium 137.3 42.86	57 La 镧 Lanthanum 138.9 46.20	58 Ce 铈 Cerium 140.1 52.09	59 Pr 镨 Praseodymium 140.9 56.48	60 Nd 钕 Neodymium 144.2 60.66	61 Pm 钷 Promethium (147) 62.94	62 Sm 钐 Samarium 150.4 67.5	63 Eu 铕 Europium 152.0 69.72	64 Gd 钆 Gadolinium 157.2 70.16	65 Tb 铽 Terbium 158.9 72.18	66 Dy 镝 Dysprosium 162.5 74.23	67 Ho 铥 Holmium 164.9 76.32	68 Er 铒 Erbium 167.2 78.46	69 Tm 铥 Thulium 168.9 80.64	70 Yb 镱 Ytterbium 173.0 82.90	71 Lu 镥 Lutetium 174.9 85.07	72 Hf 铪 Hafnium 178.4 87.90	73 Ta 钽 Tantalum 180.9 89.90	74 W 钨 Tungsten 183.8 90.94	75 Re 铼 Rhenium 186.2 93.04	76 Os 锇 Osmium 190.2 95.04	77 Ir 铱 Iridium 192.2 97.22	78 Pt 铂 Platinum 195.0 99.44	79 Au 金 Gold 197.0 101.07	80 Hg 汞 Mercury 200.5 102.90	81 Tl 铊 Thallium 204.3 104.07	82 Pb 铅 Lead 207.2 106.33	83 Bi 铋 Bismuth 208.9 108.36	84 Po 钋 Polonium (209) 111.0	85 At 砹 Astatine (210) 112.75	86 Rn 氡 Radon (222) 118.97	87 Fr 钫 Francium (223) 120.91	88 Ra 镭 Radium (226) 125.08	89 Ac 锕 Actinium (227) 137.0	90 Th 钍 Thorium 232.0 153.9	91 Pa 镤 Protactinium 231.0 158.9	92 U 铀 Uranium 238.0 183.8	93 Np 镎 Neptunium (237) 200.9	94 Pu 钚 Plutonium (244) 208.9	95 Am 镅 Americium (243) 227.0	96 Cm 锔 Curium (247) 247.0	97 Bk 锫 Berkelium (247) 247.0	98 Cf 锿 Californium (251) 251.0	99 Es 镄 Einsteinium (252) 252.0	100 Fm 镆 Fermium (257) 257.0

■ 碱金属 Alkali Metals ■ 非金属 Non-Metal ■ 过渡金属 Transitional element
■ 卤族元素 Halogen ■ 镧系元素 Lanthanoids ■ 锕系元素 Actinides
■ 碱土金属 Alkaline earth ■ 主族金属 Main group metal ■ 稀有气体 Rare gases

注1: #表示放射性元素(Radioactive Elements), *表示人造元素(Man-made Elements)
 注2: 元素名称下的数字顺序依次代表: 原子量, K α , K β , L α , L β , L γ , L ϵ
 (All the Numbers are ordered one by one in this way: Atomic Number, Element Symbol, Atomic Weight, K α , K β , L α , L β , L γ , L ϵ)

Ln	57 La 镧 Lanthanum 138.9 33.30	58 Ce 铈 Cerium 140.1 34.57	59 Pr 镨 Praseodymium 140.9 35.86	60 Nd 钕 Neodymium 144.2 37.19	61 Pm 钷 Promethium (147) 38.94	62 Sm 钐 Samarium 150.4 39.91	63 Eu 铕 Europium 152.0 41.32	64 Gd 钆 Gadolinium 157.2 42.76	65 Tb 铽 Terbium 158.9 44.23	66 Dy 镝 Dysprosium 162.5 45.73	67 Ho 铥 Holmium 164.9 47.26	68 Er 铒 Erbium 167.2 48.82	69 Tm 铥 Thulium 168.9 50.41
Actinoid	89 Ac 锕 Actinium (227) 89.79	90 Th 钍 Thorium 232.0 103.3	91 Pa 镤 Protactinium 231.0 12.85	92 U 铀 Uranium 238.0 15.71	93 Np 镎 Neptunium (237) 18.41	94 Pu 钚 Plutonium (244) 11.71	95 Am 镅 Americium (243) 15.75	96 Cm 锔 Curium (247) 18.41	97 Bk 锫 Berkelium (247) 10.77	98 Cf 锿 Californium (251) 11.05	99 Es 镄 Einsteinium (252) 113.3	100 Fm 镆 Fermium (257) 116.2	101 Md 镅 Mendelevium (258) 118.9